



# Software Requirements Specification (SRS)

**Project Name:** KidDose / DoseDek (Pediatric Dosage & Safety Calculator)

**Document Version:** 1.0.0

**Status:** Draft / System Specification

**Standard:** ISO/IEC/IEEE 29148 / IEEE 830 Compliant

## 1. Introduction

### 1.1 Purpose

This Software Requirements Specification (SRS) document describes the complete functional, non-functional, and interface requirements for the **KidDose / DoseDek** mobile application. It serves as the definitive reference for developers, quality assurance engineers, medical content reviewers, and project stakeholders.

### 1.2 Scope

KidDose / DoseDek is a cross-platform (iOS and Android) offline-first mobile application designed to assist healthcare professionals (physicians, interns, pediatricians, pharmacists, emergency staff) in calculating pediatric medication dosages accurately and instantly.

**Key capabilities include:**

- Real-time dosage calculation based on a single patient weight input ( $\text{sub-100 ms}$  rendering time).
- Automated safety checks including age contraindications, disease-specific alerts (e.g., G6PD deficiency, renal impairment), and daily/single dose caps.
- Offline-first local database execution with bundled drug data.

## 1.3 Definitions, Acronyms, and Abbreviations

- **SRS:** Software Requirements Specification
- **TDD:** Technical Design Document
- **PRD:** Product Requirements Document
- **G6PD:** Glucose-6-Phosphate Dehydrogenase Deficiency
- **OPD:** Outpatient Department
- **ER:** Emergency Room
- **Max Cap:** The maximum allowable dose (single or daily) enforced to prevent toxicity.
- **UI/UX:** User Interface / User Experience

## 1.4 Document Overview


The remainder of this document is organized into:

- **Section 2:** Overall Description (product perspective, user classes, operating environment).
- **Section 3:** Functional Requirements (detailed behavior by feature module).
- **Section 4:** External Interface Requirements (UI, hardware, software).
- **Section 5:** Non-Functional Requirements (performance, safety, security, quality attributes).
- **Section 6:** Verification & Acceptance Criteria.

# 2. Overall Description

## 2.1 Product Perspective

KidDose is a standalone, client-side mobile application with zero runtime server dependency for core calculations. It replaces manual reference books, paper dosage tables, and multi-step digital calculators with a single-input, highly responsive interface.



## 2.2 Product Functions Overview

1. **Weight & Patient Context Management:** Capture weight, age, and clinical risk flags.
2. **Dynamic Dose & Volume Calculation Engine:** Calculate target mg and liquid volume (mL).
3. **Automated Safety Evaluation:** Validate age limits, contraindications, and dose caps.
4. **Preset Category Navigation:** Filter drugs by clinical context (OPD, ER, Asthma, Antibiotics).
5. **Concentration Switching:** Toggle available liquid preparations without navigating away.
6. **Local Persistence:** Retain user preferences and favorite pinned drugs locally.

## 2.3 User Classes and Characteristics

- **General Practitioners / Interns:** High workload, fast decision-making environment. Require maximum speed, zero friction, and clear safety guardrails.

- **Pediatric Specialists:** Require high precision, access to second-line drugs, and advanced dosage adjustments.
- **ER Doctors / Nurses:** Require rapid access to emergency/resuscitation drug protocols during critical resuscitation events.
- **Pharmacists:** Require verification of concentration matching and maximum daily dose safety limits.

## 2.4 Operating Environment

- **iOS:** iOS 14.0 or higher.
- **Android:** Android 8.0 (API Level 26) or higher.
- **Connectivity:** 100% operational without network connection. Optional Wi-Fi/Cellular connectivity for database updates only.

## 2.5 Design & Implementation Constraints

- **Platform:** Developed using Flutter (Dart framework).
- **Footprint:** App binary package size must not exceed 30 MB.
- **Data Privacy:** Zero storing or transmitting of Patient Health Information (PHI/PII). All calculation inputs are transient in memory.

# 3. Specific Functional Requirements

## 3.1 Module 1: Patient Context & Weight Input

- **FR-1.1:** The app **SHALL** provide a prominent single numeric input field for Patient Weight in kilograms ( $\text{kg}$ ) on the primary dashboard.
- **FR-1.2:** The weight input **SHALL** accept floating-point numbers with up to 2 decimal places (e.g.,  $8.45 \text{ kg}$ ).
- **FR-1.3:** The weight input field **SHALL** automatically gain focus or accept rapid numpad entries without requiring modal dialogs.
- **FR-1.4:** The app **SHALL** provide optional secondary input controls for Patient Age (in months or years).
- **FR-1.5:** The app **SHALL** provide quick toggle switches for clinical risk factors:

- **[G6PD Deficiency]** (Boolean: On/Off)
- **[Renal Impairment]** (Boolean: On/Off)
- **FR-1.6:** Changing any input field **SHALL** immediately update all rendered drug results in  $\leq 100\text{ ms}$ .

## 3.2 Module 2: Dosage & Volume Calculation Engine

- **FR-2.1:** For a given drug and patient weight  $W$  ( $\text{kg}$ ), the engine **SHALL** calculate the raw single dose  $D_{\text{raw}}$  ( $\text{mg}$ ) as:  

$$D_{\text{raw}} = W \times \text{Dose}_{\text{std}}$$
 where  $\text{Dose}_{\text{std}}$  is the standard dosage parameter ( $\text{mg/kg/dose}$ ).
- **FR-2.2:** The engine **SHALL** clamp the calculated dose  $D_{\text{final}}$  ( $\text{mg}$ ) against the drug's Maximum Single Dose Cap  $M_{\text{single}}$  ( $\text{mg}$ ):  

$$D_{\text{final}} = \min(D_{\text{raw}}, M_{\text{single}})$$
- **FR-2.3:** The engine **SHALL** calculate the liquid volume  $V$  ( $\text{mL}$ ) based on selected concentration parameters  $C_{\text{mg}}$  ( $\text{mg}$ ) per  $C_{\text{ml}}$  ( $\text{mL}$ ):  

$$V = \frac{D_{\text{final}}}{C_{\text{mg}}} \times C_{\text{ml}}$$
- **FR-2.4:** Liquid volume results **SHALL** be rounded to 1 decimal place (e.g.,  $3.75\text{ mL} \rightarrow 3.8\text{ mL}$ ) for clinical clarity.
- **FR-2.5:** The engine **SHALL** evaluate total daily dosage  $D_{\text{daily}}$  against Maximum Daily Dose Cap  $M_{\text{daily}}$  ( $\text{mg}$ ):  

$$D_{\text{daily}} = D_{\text{final}} \times N_{\text{doses}}$$
 If  $D_{\text{daily}} > M_{\text{daily}}$ , the app **SHALL** flag a "Daily Dose Exceeded" warning.

## 3.3 Module 3: Real-Time Safety & Red Flag Rules

- **FR-3.1 Age Limit Check:** If  $\text{Age}_{\text{patient}} < \text{Age}_{\text{min\_drug}}$ , the drug card **SHALL** display a prominent **RED ALERT** banner stating: *"Contraindicated in patients under X months/years"*.

- **FR-3.2 G6PD Check:** If `[G6PD Deficiency]` is toggled **ON** and `$Drug.g6pdSafe == \text{false}$`, the app **SHALL** visually flag the drug with a **RED WARNING** tag and disable standard volume display until acknowledged.
- **FR-3.3 Renal Check:** If `[Renal Impairment]` is toggled **ON** and `$Drug.renalSafe == \text{false}$`, the app **SHALL** display an **AMBER WARNING** banner with dose adjustment guidance.
- **FR-3.4 Max Dose Cap Alert:** When `$D_{\text{raw}} \geq M_{\text{single}}$`, the UI **SHALL** display a "**MAX DOSE CAPPED**" indicator alongside the capped value to prevent accidental over-prescription in larger children.

### 3.4 Module 4: Category Navigation & Drug Selection

- **FR-4.1:** The UI **SHALL** provide top-level category tabs for quick filtering:
  - `OPD Common` (Default)
  - `Antibiotics`
  - `ER / Emergency`
  - `Asthma / Respiratory`
  - `Favorites / Pinned`
- **FR-4.2:** Users **SHALL** be able to pin/unpin frequently used drugs to their `Favorites` tab with a single tap.
- **FR-4.3:** The app **SHALL** include a real-time instant search bar capable of filtering drugs by Generic English name (e.g., *Amoxicillin*) and Thai name (e.g., *อะม็อกซิซิลลิน*).

### 3.5 Module 5: Concentration Selection

- **FR-5.1:** Each drug card **SHALL** display its default concentration label (e.g., `120 mg/5 mL`).
- **FR-5.2:** If multiple concentrations exist for a drug, the card **SHALL** provide an inline dropdown/toggle allowing the user to switch preparations instantaneously.
- **FR-5.3:** Selecting a new concentration **SHALL** immediately recalculate the volume (`$V$`) without affecting the weight input or other drugs.

## 4. External Interface Requirements

### 4.1 User Interface (UI) Requirements

- **Color System:**
  - **Primary Action:** Medical Blue ( #1E88E5 )
  - **Safety Red Flag:** Crimson Red ( #D32F2F )
  - **Caution / Warning:** Amber Yellow ( #FFA000 )
  - **Success / Normal:** Emerald Green ( #2E7D32 )
- **Typography:** Clean sans-serif fonts (e.g., Inter / Sukhumvit Set) with high contrast for rapid legibility under dim hospital lighting.
- **Layout:** Single-screen dashboard layout with dynamic scrolling drug cards below fixed top controls.

### 4.2 Software Interfaces

- **Local Storage Interface:** Interfacing with embedded NoSQL DB (Isar / Hive) for fetching drug schemas and user preference states.
- **Asset Manager:** Local file reader interface for reading bundled `drugs_v1.json` on cold boot.

### 4.3 Hardware Interfaces

- **Supported Screen Sizes:** Handheld smartphones ( $4.7\text{ inches}$  to  $6.8\text{ inches}$ ) and medical tablets ( $8.0\text{ inches}$  to  $12.9\text{ inches}$ ).
- **Orientation:** Portrait mode prioritized for handheld single-thumb operation.

## 5. Non-Functional Requirements

### 5.1 Performance Requirements

| Metric                     | Target Requirement  | Strict Upper Limit  |
|----------------------------|---------------------|---------------------|
| Weight Input Reaction Time | $\leq 10\text{ ms}$ | $\leq 50\text{ ms}$ |

| Metric                      | Target Requirement          | Strict Upper Limit       |
|-----------------------------|-----------------------------|--------------------------|
| <b>Cold App Launch Time</b> | $\leq 800 \text{ ms}$       | $\leq 1200 \text{ ms}$   |
| <b>Warm App Launch Time</b> | $\leq 200 \text{ ms}$       | $\leq 400 \text{ ms}$    |
| <b>Scroll Frame Rate</b>    | $60 \text{ fps}$ continuous | $50 \text{ fps}$ minimum |
| <b>Memory Consumption</b>   | $\approx 35 \text{ MB}$     | $\leq 60 \text{ MB}$     |

## 5.2 Safety & Reliability Requirements

- **Calculation Precision:** Floating-point arithmetic must maintain accuracy to 4 decimal places internally before rendering.
- **Zero Failure Rate on Red Flags:** All safety flag conditions (\$Age\$, \$G6PD\$, \$Renal\$, \$MaxCap\$) must execute 100% deterministically before UI rendering occurs.
- **Data Validation:** Negative weight inputs or non-numeric characters must be blocked at input level.

## 5.3 Security & Compliance Requirements

- **Zero PHI Storage:** The system shall not log, store, or transmit patient names, weight history, or clinical choices.
- **Offline Independence:** No external network requests allowed during dosage calculation to ensure uninterrupted operation during network outages.

## 5.4 Usability & Accessibility Requirements

- **Click Count Target:** Maximum of  $1 \text{ tap}$  required (entering weight) to view dosages for standard OPD drugs.
- **Touch Targets:** Minimum touch target size of  $48 \times 48 \text{ dp}$  for all interactive buttons and toggles.

# 6. Verification & Traceability Matrix



| Requirement ID | Description                               | Verification Method    | Acceptance Criteria                                                                                                        |
|----------------|-------------------------------------------|------------------------|----------------------------------------------------------------------------------------------------------------------------|
| FR-1.1         | Single Weight Input Field                 | Inspection / Test      | Field visible on cold launch; accepts numeric values.                                                                      |
| FR-1.6         | Real-time Update<br>$\leq 100 \text{ ms}$ | Automated UI Benchmark | UI recalculates and renders all cards in $< 50 \text{ ms}$ upon typing.                                                    |
| FR-2.2         | Max Single Dose Cap                       | Unit Test              | Patient $60 \text{ kg}$ on Paracetamol ( $10 \text{ mg/kg}$ ) yields capped $500 \text{ mg}$ instead of $600 \text{ mg}$ . |
| FR-3.1         | Age Restriction Flag                      | Unit Test & UI Test    | Setting age to 14 months for Chlorpheniramine displays Red Flag warning.                                                   |
| FR-3.2         | G6PD Warning Trigger                      | Unit Test & UI Test    | Toggling G6PD ON displays red warning on unsafe antimalarial / antibiotic cards.                                           |
| NFR-5.3        | Offline Execution                         | Integration Test       | All functions remain fully operational in Airplane Mode.                                                                   |

## 7. Medical Disclaimer Specification

The application **MUST** display a non-intrusive legal disclaimer during first launch:

**Medical Disclaimer:** *KidDose / DoseDek is designed as a clinical decision-support tool for licensed healthcare professionals. It does not replace clinical judgment. Prescribers are solely responsible for verifying final dosages against institutional protocols and patient specifics.*